

November 2010

Lichtleitkabel
Fibre lightguide

MEDI-LIGHT HTX

Hoch
Temperaturbeständig
Xenon geeignet



1,8 / 2,3 / 3,0 / 3,5 m

NEU / NEW

Spezialfiberglas

- hohe Lichtübertragung
- stets weißes Licht

Knickfest

- lange Lebensdauer
bei optimalem Licht

Flexibel

- für große Beweglichkeit
und geringe Ermüdung

Aufbereitung und Sterilisation

- waschmaschinenbeständig
und autoklavierbar bei 134 °C

Universell

- adaptierbar an alle gebräuchlichen
Endoskope und Projektoren,
auch XENON

Special fiberglass

- intense light transfer
- constant white light

Resistant to buckling

- durable with optimal light

Very light and flexible

- allowing great movability
without fatigue

Cleaning and sterilization

- is suitable for washing machines
and can be sterilized
in an autoclave 134 °C

Universal

- can be adapted to all usual
endoscopes and projectors

Lichtleitkabel Fibre lightguide

MEDI-LIGHT HTX

Lichtleitkabel / fibre lightguide

Kabel allein, Adapterauswahl siehe unten / Cable only, choose adapter



Ø 3,5 mm

Ø 4,8 mm

Hoch
Temperaturbeständig
Xenon geeignet

Länge / length

Art.-No.

Art.-No.

1,8 m

85.056.18

85.058.18

2,3 m

85.056.23

85.058.23

3,0 m

85.056.30

85.058.30

3,5 m

85.058.35

Adapter / Adapters



Lichtquelle
Lightsource



Endoskop
Instrument


Adapter passend für

Adapter passend für



85.060.01
- Storz / Martin /
Wisap / Aesculap

85.060.51
- Storz
- Olympus
- Winter & Ibe
- Aesculap
- Martin
- Wisap



85.060.02
- Wolf

85.060.52
- Wolf



85.060.03
- Olympus

... weitere Lichtleitkabel und Adapter nach Ihren speziellen Wünschen
... other lightguides and adapters according to your special request