

Januar 2013

Fibre lightguide

MEDI-LIGHT HTX

Hoch
high
Temperaturbeständig
temperature resistant
Xenon geeignet
suitable for Xenon

1,8 / 2,3 / 3,0 / 3,5 m



Special fiberglass

- intense light transfer
- constant white light

Resistant to buckling

- durable with optimal light

Very light and flexible

- allowing great movability
without fatigue

Cleaning and sterilization

- is suitable for washing machines
and can be sterilized
in an autoclave 134 °C

Universal

- can be adapted to all usual
endoscopes and projectors

Fibre lightguide MEDI-LIGHT HTX

Fibre lightguide

Cable only, choose adapter



Hoch
high
Temperaturbeständig
temperatur resistant
Xenon geeignet
suitable for Xenon

Length

- 1,8 m
- 2,3 m
- 3,0 m
- 3,5 m

Ø 3,5 mm

Art.-No.

- 85.056.18
- 85.056.23
- 85.056.30
-

Ø 4,8 mm

Art.-No.

- 85.058.18
- 85.058.23
- 85.058.30
- 85.058.35

Adapters



Lightsource



Instrument

Adapters suitable for



- 85.060.01**
- Storz / Martin /
Wisap / Aesculap




- 85.060.11**
- Storz



- 85.060.02**
- Wolf




- 85.060.03**
- Olympus



- 85.060.07**
- Fuji 120-er Serie

Adapters suitable for

- 85.060.51**
- Storz
- Olympus
- Winter & Ibe
- Aesculap
- Martin
- Wisap



- 85.060.52**
- Wolf



- 85.060.53**
- ACM-British Standard
- Heine
- Welch Allyn



... other lightguides and adapters according to your special request