

Juni 2022

Bipolar forceps

compatible to:

- Berchtold / Integra
- BOWA
- Erbe
- Martin
- Valleylab



		
	needle	0,5 mm
Length	pointed	very fine

straight



11 cm	90.822.11 (H-HF9611-11)	90.824.11 (H-HF9612-11)
-------	-----------------------------------	-----------------------------------

angled



11 cm	90.825.11 (H-HF9621-11)	90.827.11 (H-HF9622-11)
-------	-----------------------------------	-----------------------------------

Bipolar forceps




			
Length	0,5 mm pointed	1 mm blunt	2 mm large

straight

16 cm	90.822.16 (H-HF9612-16)	90.824.16 (H-HF9613-16)	90.826.16 (H-HF9614-16)
19,5 cm		90.824.20 (H-HF9613-19)	90.826.20 (H-HF9614-19)
22 cm		90.824.22 (H-HF9613-22)	90.826.22 (H-HF9614-22)

Bipolar forceps




			
Length	0,5 mm pointed	1 mm blunt	2 mm large

angled

16 cm	90.825.16 (H-HF9622-16)	90.827.16 (H-HF9623-16)	90.829.16 (H-HF9624-16)
19,5 cm		90.827.20 (H-HF9623-19)	90.829.20 (H-HF9624-19)
22 cm		90.827.22 (H-HF9623-22)	90.829.22 (H-HF9624-22)

Bipolar forceps




			
Length	0,5 mm pointed	1 mm blunt	2 mm large

bayonet

16 cm	90.842.16 (H-HF9632-16)	90.844.16 (H-HF9633-16)	
19,5 cm		90.844.20 (H-HF9633-19)	90.846.20 (H-HF9634-19)
22 cm		90.844.22 (H-HF9633-22)	90.846.22 (H-HF9634-22)
24 cm		90.844.24 (H-HF9633-24)	90.846.24 (H-HF9634-24)


Bipolar forceps



bayonet, angled up



19,5 cm	90.845.20 (H-HF9653-19)
22 cm	90.845.22 (H-HF9643-22)



bayonet, angled down

19,5 cm	90.851.20 (H-HF9643-19)
---------	-----------------------------------


Bipolar forceps



Length	1,0 mm blunt

curved

19 cm	90.859.19 (H-HF9693-19)
-------	-----------------------------------



max
500
Vp











Length	1 mm blunt

approx. 60° curved

17 cm	90.859.18 (H-HF9693-17)
-------	-----------------------------------

Cables for bipolar forceps



Electrosurgical unit		Art.-No. / Length	Instrument
Erbe ACC / ICC T-Modell		90.286.58 = 4 m	 "P-Type" angled
Berchtold/Integra Martin Aesculap		90.287.58 = 4 m	
Valleylab		90.287.38 = 5 m	
Erbe ACC / ICC T-Modell		90.286.56 = 4 m	 "P-Type"
Berchtold/Integra Martin Aesculap		90.287.56 = 4 m	
Erbe ACC / ICC T-Modell		90.286.30 = 3 m 90.286.50 = 5 m	 Standard
Berchtold/Integra Martin Aesculap		90.287.30 = 3 m 90.287.50 = 5 m	
Valleylab		90.287.37 = 5 m	

More cables see leaflet *051126-88



= Packaging unit



= Max. temperature
in autoclave



= max. rated
voltage

Numerical Index

Art.-No./	Page	Art.-No./	Page
90.286.30	7	90.829.16	3
90.286.50	7	90.829.20	3
90.286.56	7	90.829.22	3
90.286.58	7	90.842.16	4
90.287.30	7	90.844.16	4
90.287.37	7	90.844.20	4
90.287.38	7	90.844.22	4
90.287.50	7	90.844.24	4
90.287.56	7	90.845.20	5
90.287.58	7	90.845.22	5
90.822.11	1	90.846.20	4
90.822.16	2	90.846.22	4
90.824.11	1	90.846.24	4
90.824.16	2	90.851.20	5
90.824.20	2	90.859.18	6
90.824.22	2	90.859.19	6
90.825.11	1		
90.825.16	3		
90.826.16	2		
90.826.20	2		
90.826.22	2		
90.827.11	1		
90.827.16	3		
90.827.20	3		
90.827.22	3		