

Mai 2017

## Autopsy saw standard type



**30.218.20** Oscillating autopsy saw, 230 V~ 50/60 Hz and 500 Watt.

- 1 piece **30.210.50** sawblade Ø50mm
- 1 piece **30.210.65** sawblade Ø65mm
- 1 piece **30.212.09** sawblade for autopsy, bilateral
- 1 pair of flat wrenches **30.212.02** for the exchange of the sawblades

### Technical data:


**Item number:** **30.218.20**

Connection for power supply: 230 V~ 50/60 Hz  
Input: 500 Watt  
Adjustment range: 3500 to 24000 osz./min. load independent


Protection class: II according to VDE 0740 and respectively BF  
Weight: 1.500gr.  
Mains cable length: 5m, EU-connection

According to EMV and EN60601-1-1, EN60601-1-2, IEC 60601-1-1, IEC 60601-1-2

## Sawblades for autopsy saws



**30.210.50**  
Sawblade Ø50mm




**30.210.65**  
Sawblade Ø65mm




**30.212.09**  
Sawblade bilateral



**30.219.84**  
Depth stop eccentric



**30.212.05**  
Segmented saw blade Ø65mm,  
40mm of width (Butterfly)



**30.212.08**  
Depth saw blade 20mm of width

**30.212.06**  
Depth saw blade 8mm of width