

Mai 2017

Autopsy saw standard type



30.218.20 Oscillating autopsy saw, 230 V~ 50/60 Hz and 500 Watt.

- 1 piece **30.210.50** sawblade Ø50mm
- 1 piece **30.210.65** sawblade Ø65mm
- 1 piece **30.212.09** sawblade for autopsy, bilateral
- 1 pair of flat wrenches **30.212.02** for the exchange of the sawblades

Technical data:

Item number: **30.218.20**

Connection for power supply: 230 V~ 50/60 Hz

Input: 500 Watt

Adjustment range: 3500 to 24000 osz./min. load independent


Protection class: II according to VDE 0740 and respectively BF

Weight: 1.500gr.


Mains cable length: 5m, EU-connection

According to EMV and EN60601-1-1, EN60601-1-2, IEC 60601-1-1, IEC 60601-1-2

Sawblades for autopsy saws



30.210.50
Sawblade Ø50mm




30.210.65
Sawblade Ø65mm




30.212.09
Sawblade bilateral



30.219.84
Depth stop eccentric



30.212.05
Segmented saw blade Ø65mm,
40mm of width (Butterfly)



30.212.08
Depth saw blade 20mm of width

30.212.06
Depth saw blade 8mm of width