

Juli 2017

**Autopsy saw for use with triangular shafts
with detachable saw head for disinfection and autoclaving**



30.219.20 Oscillating autopsy saw, 230 V~ 50/60 Hz and 500 Watt. The saw head is detachable and can be disinfected and autoclaved.

Standard equipment:

- 1 piece **30.219.92** shaft 39mm (used length 5mm) for saw blades
- 1 piece **30.210.50** saw blade Ø50mm
- 1 piece **30.210.65** saw blade Ø65mm
- 1 piece **30.212.09** segmented saw blade Ø65mm, 39mm of length, bilateral
- 1 piece **30.212.10** T-handle to detach the oscillating head and triangular shafts
- 1 pair **30.219.96** wrenches SW 14mm for the exchange of the saw blades

Technical data:

Item number:	30.219.20
Connection for power supply:	230 V~ 50/60 Hz
Input:	500 Watt
Adjustment range:	3500 to 24000 osz./min. load independent
Protection class:	II according to VDE 0740 and respectively BF
Weight:	1.500gr.
Mains cable length:	5m, EU-connection

According to EMV and EN60601 ff.


**Autopsy saw for use with triangular shafts
with detachable saw head for disinfection and autoclaving**

**Stainless saw blades for
autopsy saws**

Single / Replacement parts



30.219.80
Double saw Ø45mm consisting of:
shaft 30.219.91, saw blade 30.219.85
and adjustable saw blade 30.219.81



30.219.85
Saw blade Ø45mm for double saw 30.219.80



30.219.81
Saw blade Ø45mm adjustable
for double saw 30.219.80



30.219.70
Saw blade Ø65mm bilateral for autopsy
consisting of: shaft 39mm 30.219.92 and
bilateral saw blade 30.212.09




30.212.09
Saw blade bilateral

**Autopsy saw for use with triangular shafts
with detachable saw head for disinfection and autoclaving**



30.219.50

Saw blade Ø50mm with shaft consisting of: shaft 39mm 30.219.92 and saw blade 30.210.50



30.210.50

Saw blade Ø50mm

30.210.65

Saw blade Ø65mm

30.219.65

Saw blade Ø65mm with shaft consisting of: shaft 39mm 30.219.92 and saw blade 30.210.65



30.219.83

Saw blade Ø65mm with depth stop consisting of: shaft 58mm 30.219.90, saw blade 30.210.65 and depth stop 30.219.84




30.219.84

Depth stop eccentric

**Autopsy saw for use with triangular shafts
with detachable saw head for disinfection and autoclaving**




30.219.71
Segmented saw blade Ø65mm, 40mm
of width (Butterfly) with shaft consisting of:
shaft 39mm 30.219.92 and saw blade
30.212.05



30.212.05
Segmented saw blade Ø65mm,
40mm of width (Butterfly)



30.219.75
Depth saw blade 20mm of width with shaft
consisting of: shaft 39mm 30.219.92
and saw blade 30.212.08



30.212.08
Depth saw blade 20mm of width

30.219.76
Depth saw blade 8mm of width with shaft
consisting of: shaft 39mm 30.219.92
and saw blade 30.212.07

30.212.06
Depth saw blade 8mm of width



30.219.90
Shaft 58mm
(used length 30mm)



30.219.92
Shaft 39mm
(used length 5mm)



30.219.91
Shaft for
double saw blade