


September 2018

Bühler-*ErgonoMIC*

Range of instruments for endoscopy

5 mm

- for outstanding precision during endoscopic procedures
- axial handle permits easy rotation, e.g. for endoscopic needle holders with tungsten carbide inserts
- safe use due to durable construction and reliable mechanics
- easy to clean and maintain due to simple disassembly
- 5 year guarantee * - cost-effective due to long life for instrument




* on defects in material and manufacturing

Page 1


Bühler-ErgonoMIC

Range of instruments for endoscopy


5 mm



Outstanding precision and reliable mechanics for highest performance during endoscopic procedures



Precise control due to easy rotation, e.g. during suturing with a needle holder equipped with tungsten carbide inserts



Easy to clean and maintain due to simple disassembly

Bühler-*ErgonoMIC*

Range of instruments for endoscopy

5 mm



Original size

M: 1:1


85.405.05

Bühler-*ErgonoMIC* TC needle holder, 5mm, straight jaws, single action, jaws with tungsten carbide insert for reliable grasping of the needle

Bühler-*ErgonoMIC*

Range of instruments for endoscopy

5 mm




Original size

M: 1:1

85.405.15

Bühler-*ErgonoMIC* "Small" grasper, 5mm, single action, a universal grasping forceps




Original size

M: 1:1

85.405.16

Bühler-*ErgonoMIC* "ALLIS/KOCHER" grasper, 5mm, single action, for grasping tissues with the tip of the jaws, esp. for inflamed organs



Original size

M: 1:1


85.405.23

Bühler-*ErgonoMIC* "ALLIGATOR" grasper, 5mm, single action, aggressive grasper, also known as spiked jaws

Bühler-*ErgonoMIC*

Range of instruments for endoscopy

5 mm




Original size

M: 1:1

85.405.25

Bühler-*ErgonoMIC* “Spoon-shaped/biopsy” grasper, 5mm, single action, for biopsy or to grasp small stones




Original size

M: 1:1

85.405.29

Bühler-*ErgonoMIC* “Atrauma” grasper, 5mm, single action, for gently grasping delicate structures thru a non-traumatic De Bakey serration



Original size

M: 1:1

85.405.45

Bühler-*ErgonoMIC* “Non-traumatic intestinal” grasper, 5mm, single action, for gently grasping of larger tissue structures thru a non-traumatic De Bakey serration

Bühler-*ErgonoMIC*

Range of instruments for endoscopy

5 mm

Notes

Bühler-ErgonoMIC


Range of instruments for endoscopy

5 mm



85.019.85

Positioning rack for M.I.S instruments with diameters of 5mm and 10mm, autoclavable, rack can be placed into wire mesh trays




***960529-33**

Silicon mat 540x240mm, blue




88.691.05

Sterilizing tray, wire mesh 540x255x50mm



85.405.97 = 5mm

Outer sleeve for ErgonoMIC instruments



85.405.98 = 5mm

Coupling nut for ErgonoMIC instruments, golden

Bühler-ErgonoMIC

Range of instruments for endoscopy

5 mm

Numerical Index

Cat.-No.	Page
85.019.85	7
85.405.05	3
85.405.15	4
85.405.16	4
85.405.23	4
85.405.25	5
85.405.29	5
85.405.45	5
85.405.97	7
85.405.98	7
85.691.05	7
*960529-33	7