

April 2020

Finger ring saw with diamond saw blade

Our finger ring saw with diamond-tipped saw blade allows finger rings made of gold and silver to be cut through easily and quickly. Equipped with a specially coated (diamond) saw blade, it is **also ideal for finger rings made of stainless steel.**



Item No. *100706-33

The ergonomically designed handle ensures the best power transmission to the diamond saw blade.

The instrument is made of highly hardened surgical stainless steel.

The wide wings of the rotary saw blade handle also enable high power transmission.

The direction of rotation of the saw blade is clearly indicated with a large arrow.

An ideal tool for finger ring removal not only for the emergency room, but also indispensable for every jeweler.

April 2020

Finger ring saw with diamond saw blade

Saw disassembled:



Close-up of the diamond saw blade:



: