

Januar 2013

*Fibre lightguide*

**MEDI-LIGHT HTX**

**Hoch**  
high  
**Temperaturbeständig**  
temperature resistant  
**Xenon geeignet**  
suitable for Xenon

1,8 / 2,3 / 3,0 / 3,5 m



**Special fiberglass**

- intense light transfer
- constant white light

**Resistant to buckling**

- durable with optimal light

**Very light and flexible**

- allowing great movability  
without fatigue

**Cleaning and sterilization**

- is suitable for washing machines  
and can be sterilized  
in an autoclave 134 °C

**Universal**

- can be adapted to all usual  
endoscopes and projectors

**Fibre lightguide MEDI-LIGHT HTX**

**Fibre lightguide**

*Cable only, choose adapter*



Hoch  
high  
Temperaturbeständig  
temperatur resistant  
Xenon geeignet  
suitable for Xenon

**Length**

- 1,8 m
- 2,3 m
- 3,0 m
- 3,5 m

**Ø 3,5 mm**

**Art.-No.**

- 85.056.18
- 85.056.23
- 85.056.30
- 

**Ø 4,8 mm**

**Art.-No.**

- 85.058.18
- 85.058.23
- 85.058.30
- 85.058.35

**Adapters**



*Lightsource*

*Instrument*




*Adapters suitable for*

*Adapters suitable for*



**85.060.01**  
- Storz / Martin /  
Wisap / Aesculap

- 85.060.51**  
- Storz  
- Olympus  
- Winter & Ibe  
- Aesculap  
- Martin  
- Wisap



**85.060.11**  
- Storz




**85.060.02**  
- Wolf

**85.060.52**  
- Wolf



**85.060.03**  
- Olympus

- 85.060.53**  
- ACM-British Standard  
- Heine  
- Welch Allyn



**85.060.07**  
- Fuji 120-er Serie

*... other lightguides and adapters according to your special request*