

April 2019

Fibre lightguide

MEDI-LIGHT HTX

Hoch
high
Temperaturbeständig
temperature resistant
Xenon geeignet
suitable for Xenon

1,8 / 2,3 / 3,0 / 3,5 m



Special fiberglass

- intense light transfer
- constant white light

Resistant to buckling

- durable with optimal light

Very light and flexible

- allowing great movability
without fatigue

Cleaning and sterilization

- is suitable for washing machines
and can be sterilized
in an autoclave 134 °C

Universal

- can be adapted to all usual
endoscopes and projectors

April 2019

Fibre lightguide **MEDI-LIGHT HTX**

Fibre lightguide

Cable only, choose adapter



Hoch
high
Temperaturbeständig
temperatur resistant
Xenon geeignet
suitable for Xenon

Length

1,8 m
2,3 m
3,0 m
3,5 m

Ø 3,5 mm

Art.-No.

85.056.18
85.056.23
85.056.30

Ø 4,8 mm

Art.-No.

85.058.18
85.058.23
85.058.30
85.058.35

Adapters



Lightsource



Instrument

Adapters suitable for

Adapters suitable for




85.060.01

- Storz
- Martin
- Wisap
- Aesculap
- Schölly

85.060.51

- Storz
- Olympus
- Winter & Ibe
- Aesculap
- Martin
- Wisap



85.060.11

- Storz

85.060.52

- Wolf



85.060.02

- Wolf

85.060.53

- ACM-British Standard
- Heine
- Welch Allyn



85.060.03

- Olympus

... other lightguides and adapters according to your special request