

Juni 2019

Fibre lightguide

MEDI-LIGHT HTX

Hoch
high
Temperaturbeständig
temperature resistant
Xenon geeignet
suitable for Xenon

1,8 / 2,3 / 3,0 / 3,5 m



Special fiberglass

- intense light transfer
- constant white light

Resistant to buckling

- durable with optimal light

Very light and flexible

- allowing great movability
without fatigue

Cleaning and sterilization

- is suitable for washing machines
and can be sterilized
in an autoclave 134 °C

Universal

- can be adapted to all usual
endoscopes and projectors

Juni 2019

***Fibre lightguide* MEDI-LIGHT HTX**

Fibre lightguide










Cable only, choose adapter



Hoch
high
Temperaturbeständig
temperatur resistant
Xenon geeignet
suitable for Xenon

<i>Length</i>	Ø 3,5 mm Art.-No.	Ø 4,8 mm Art.-No.
1,8 m	85.056.18	85.058.18
2,3 m	85.056.23	85.058.23
3,0 m	85.056.30	85.058.30
3,5 m	---	85.058.35

Adapters

 <i>Lightsource</i>	<i>Instrument</i> 
<i>Adapters suitable for</i>	<i>Adapters suitable for</i>
 <p>85.060.01 - Storz - Martin - Wisap - Aesculap - Schölly</p>	<p>85.060.51 - Storz - Olympus - Winter & Ibe - Aesculap - Martin - Wisap</p> 
 <p>85.060.11 - Storz</p>	<p>85.060.52 - Wolf</p> 
 <p>85.060.02 - Wolf</p>	<p>85.060.53 - Heine - Welch Allyn - ACM-British Standard</p> 
 <p>85.060.03 - Olympus</p>	

... other lightguides and adapters according to your special request