

Juni 2013

TEST-FOX UNIVERSAL 2.0

Testing device for electrosurgical accessories

– safe – universal – handy –



for function test of cables, electrosurgical pens and electrode pens/handles with cords etc. for popular accessories of brands:

for **Bühler**
for **Erbe**
for **Martin**

for **Berchtold**
for **Integra**
for **Aesculap**

for **Storz**
for **Valleylab**
etc.

TEST-FOX UNIVERSAL 2.0

Testing device for electrosurgical accessories

- Safety for CSSD and operating theatre due to tested RF-accessories
- Safe and easy testing without adapters
- Sophisticated and user friendly conception
- Simple handling thru a coloured identification of the testing areas
- Reasoned and evident test steps
- Accessories for easy insulation testing optionally for example for insulated lap-/chole-instruments



Testing area for monopolar cords, electrode pens/handles



Testing area for bipolar connection cables for bipolar forceps and lap-/chole-instruments



Testing area for cables for neutral electrodes



Testing area for isolations

Art.no. 91.801.00

TEST-FOX UNIVERSAL 2.0 complete with wall power supply and contacting pins Ø 2,4 mm and Ø 4,0 mm